

Is The Vial Centrifuge NRTL rated - FAQ

Date: 15-NOVEMBER-2023 Last Updated: 18-JUNE-2025

At this time, the Vial Centrifuge™ is not certified by a Nationally Recognized Testing Laboratory (NRTL) such as UL, CSA, or ETL. Therefore, it does not carry an NRTL mark for compliance with North American electrical safety standards.

However, the Vial Centrifuge™ is CE marked, indicating that it meets the essential health, safety, and environmental protection requirements outlined in applicable European Union directives. This includes compliance with standards related to:

- Low Voltage Directive (LVD)
- Electromagnetic Compatibility (EMC)
- Machinery Safety Requirements

The CE mark ensures that the product has undergone a conformity assessment and is suitable for use in regulated laboratory environments within the EU and other CE-recognizing regions.

Additional Notes for Advanced Buyers

- **Use in Regulated Facilities:** If your facility requires NRTL certification for all powered equipment (e.g., certain U.S. laboratories, universities, or government institutions), please consult your safety officer or compliance team before purchase.
- **Electrical Safety:** The Vial Centrifuge™ is designed with built-in safety features and undergoes rigorous internal quality control testing, but it has not been third-party tested under NRTL protocols.
- **Documentation:** CE Declaration of Conformity and technical specifications are available upon request for audit or procurement documentation.

👉 Click [HERE](#) for Vial Centrifuge™ ordering information and product details.