

## Is there a Difference in the Properties of Different Colors of Rubber in Autosampler Vial Cap Septa - HPLC Primer

*Date: 27-DECEMBER-2023 Last Updated: 27-OCTOBER-2025*

Different colors do not always indicate a difference in the septa properties of vial cap septa.

Septa can have slightly different properties regardless of color but the differences will always be indicated in the product Description or the Additional Information on the product page.

[Click HERE for Autosampler Vial Cap Ordering Information.](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)