



Luer Fittings Used Fr a Flash Chromatography Column - Tech Information

Date: 10-JANUARY-2024 Last Updated: 1-DECEMBER-2025

To connect the Flash Column with a female Luer inlet, male Luer outlet, to an HPLC system with a typical 1/16" SS Tubing with 10-32 threaded fittings, you will need the following:

1. [47451-42](#) Adapter: 10-32 female to male luer connection (for inlet)
2. 53.9.200124 Adapter: PP, female Luer to 10-32 male (for outlet)

ARE-APPLIED RESEARCH™

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com