

## Ferrules Choice When Using Soft Polymeric Tubing - Tips & Suggestions

*Date: 25-FEBRUARY- 2024   Last Updated: 1-DECEMBER-2025*

Kel-F™ (*PCTFE which is similar to PTFE*) and Tefzel™ (*ETFE*) ferrules give the best results when using PTFE, FEP, PFA or ETFE tubing.

[More information about PEEK tubing](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)