

Are the Chemical Solutions qualification solutions for HPLC NIST traceable - FAQ

Date: 14-JUNE-2024 Last Updated: 5JULY-2025

Yes, all Chemical Solutions™ brand PQ (Performance Qualification) and HSQ (High Speed Qualification) solutions are NIST traceable.


What Does NIST Traceable Mean?

NIST traceability ensures that the measurements and standards used in the preparation of our solutions are linked to the **National Institute of Standards and Technology (NIST)** through an unbroken chain of comparisons, all with stated uncertainties.

NIST Traceability by Product Type:

- All Chemical Solutions™ products are prepared using NIST-traceable instruments and protocols in a controlled laboratory environment.
- Our Holmium Oxide solution is directly compared to NIST specifications during manufacturing, ensuring full traceability.
- While there is no direct NIST Standard Reference Material (SRM) for our Caffeine solution, it is produced using calibrated, NIST-traceable equipment and validated procedures to maintain traceability and consistency.

These practices ensure that our solutions meet the rigorous standards required for HPLC system qualification and calibration in regulated environments.

 Click [HERE](#) for Qualification Kits Ordering Information

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com