

What is the Maximum Temperature for the Solid Bottle Caps for Mobile Phase Bottles - FAQ

Date: 31-January-2025 Last Update: 8-NOVEMBER-2025

Our mobile phase caps are made of polypropylene, and the recommended maximum temperature is 140°C.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com