



What is the Temperature Range for Non-Coated Mobile Phase Reservoir Bottles - FAQ

Date: 12-MARCH-2025 Last Updated: 8-NOVEMBER-2025

The temperature range for [non-coated Mobile Phase Reservoir Bottles](#) is -70° to 500°C. Please remember to verify the cap's temperature range as well and avoid fully sealing the bottles when heating.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com