

## Is the Zero or Controlled Flow Capillary Available in Longer than 1 Meter Length - FAQ

Date: 7-April-2025 Last Updated: 3-NOVEMBER-2025

At this time, Zero Flow and Controlled Flow capillaries are only available in 1-meter lengths.

### Why?

These specialized capillaries require a precise and uniform inner wall treatment to achieve their unique flow characteristics. Extending the length beyond 1 meter introduces significant challenges in:

- Maintaining consistent surface bonding throughout the entire capillary
- Ensuring reproducible electroosmotic flow (EOF) control
- Preserving performance and quality standards

For these reasons, longer lengths are not currently offered, as they would compromise the reliability and consistency of the product.

---

### Need a Custom Solution?

If your application requires a longer capillary or a specific configuration, please contact us. Our technical team is happy to explore potential alternatives or custom options that may suit your workflow.



Click [HERE](#) for ordering information and pictures