



What is the durability of the SOLVFIL when used with standard mobile phases and buffered mobile phases - FAQ

The longevity of filter material is influenced by the nature and quantity of substances undergoing filtration. Due to the extensive diversity of solvents—including reverse-phase solvents, normal-phase solvents, and various buffers with differing concentrations—establishing a standardized lifespan for these filters would be difficult given the multitude of variables involved. Therefore, it is advisable to replace filters when their performance noticeably decreases, such as when the filtration process becomes slower than anticipated.

Click [HERE](#) for SOLVFi™ mobile phase bottle top filter degasser ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com