



Temperature Range of MICROSOLV Brand Autosampler Vial Screw Caps - Tech Information

Date: 16-MAY-2018 Last Updated: 27-OCTOBER-2025

Most polypropylene autosampler vial caps, including MICROSOLV brand caps are designed for use at room temperature with a typical temperature range of 0°C to 60°C. It is best to stay within the specifications of the products for best results and safety.

Caution: An extreme temperature range of -55°C to 220°C is possible, however the silicone / PTFE septa will withstand a higher temperature but structural integrity of the polypropylene cap is lost.

AUTOSAMPLER VIALS AND CAPS

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com