



Are RSA Autosampler Vials Pyrogen Free - FAQ

Date: 1-JUNE-2012 Last Updated: 8-NOVEMBER-2025

YES. The RSA™ or Reduced Surface Activity autosampler vials are pyrogen free when manufactured, as they are produced sterile and are never touched by "human hands". They are then packaged, including shrink wrapped in a sterile, class 10,000 clean room environment.

We do not classify these vials as sterile for sale since the packaging is not sterilized before use. However, this can give you a very good idea of the cleanliness if you wished to sterilize them.

The entire manufacturing process is ISO Certified.

[More information about RSA vials, inserts and caps](#)



Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com