



Do you have a tear off crimp cap that is also magnetic - FAQ

No, we do not offer a magnetic tear off cap at this time because it requires aluminum (a soft metal) to allow easy tearing of the cap and aluminum is not magnetic. The magnetic material used in crimp caps is not soft enough to tear.

Click [HERE](#) for crimp cap ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com