



RSA-Pro X Surface Treatment Process Does Not Add Chlorides to the Glass Surface - Tech Information

Date: 19-MAY-2017 Last Updated: 13-DECEMBER-2025

No, the RSA-Pro X™ surface treatments used do not add chlorides to the surface of the glass.

Click [HERE](#) for the RSA-Pro X technical information page with a demonstration video.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com