



## What are the temperature specifications of the autosampler vial storage box and cap for 4ml vials - FAQ

Catalog number 9504SB-03 is temperature resistant in the range -80°C to +100°C.

Click [HERE](#) for storage box for 15x45mm autosampler vial ordering information.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)