

## What is the Purpose of a Partial Slit Septa in Autosampler Vial Caps - FAQ

*Date: 30-MARCH-2012 Last Updated: 27-OCTOBER-2025*

The partial slit septa offers both a penetration aid for thick and dull needles plus prevents evaporation before injection resulting in concentration changes. Only the silicone rubber is cut, the PTFE is not pre cut acting as an evaporation barrier.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)