



## Can I use a Vial like 95025-VN-05C with Crimp Cap like 95025-08-1S in an Oven at 105C for Approximately 1 Month - FAQ

*Date: 5-SEPTEMBER-2017 Last Update: 27-OCTOBER-2025*

**Yes, the products should performing well at this temperature and time period.**

[95025-VN-05C | MICROSOVLV Technology Corp](#)

[95025-08-1S | Crimp Cap | Silicone/PTFE | 20mm | MICROSOVLV Technology Corp](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)