

What is the Advantage of Using a Pharmafix Type of Septa for Headspace Analysis - FAQ

Date: 22-APRIL-2012 Last Updated: 27-OCTOBER-2025

The advantage is that you get a better seal with the “Pharmafix” seal. However, you should be aware that it may not be as inert as a silicone/PTFE septa designed for headspace analysis.

Click [HERE](#) for ordering information and pictures of headspace caps.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com