

## What is the Purpose of foil as Part of a Septa in Headspace Analysis - FAQ

*Date: 22-APRIL-2012 Last Updated: 27-OCTOBER-2025*

The purpose of the foil is generally to be a barrier film. To keep any solvents or gases from affecting the septa in the crimp cap. It also has a heat insulator property.

Click [HERE](#) for more information about Headspace Caps

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)