



Which Crimp Cap is Most Commonly Used for Residual Solvent Testing by the Pharmaceutical Industry - FAQ

Date: 22-APRIL-2012 Last Updated: 27-OCTOBER-2025

The crimp cap with silicone rubber / PTFE septa is most used by the pharmaceutical industry to be in compliance with USP <467>.

Our catalog number [95025-09-1S](#) is one of many caps used.

Click [HERE](#) for more information about Headspace vials and caps

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com