



Will all MICROSOLV Headspace Crimp Caps Work for all Brands of 20mm Crimp Top Vials - FAQ

Date: 22-APRIL-2012 Last Updated: 27-OCTOBER-2025

Yes, most definitely.

However, the crimp top (serum finish) for MICROSOLV vials are tested for fit and form to insure precise crimping each and every time. Always be sure to have your crimping tool calibrated approximately every 50 crimps. We cannot predict or control vials that we do not produce.

For more information about calibrating your crimping tool, please contact Technical Service for assistance.

[Click here for more information about Headspace Vials and Caps](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com