

Can BASIK Screw Caps be Used in Ion Chromatography - FAQ

Date: 29-SEPTEMBER-2014 Last Updated: 8-NOVEMBER-2025

While BASIK™ brand screw caps are not designed for ion chromatography, they may work for some applications. But I would not use them sensitive methods.

We do suggest it might be worthwhile for maximum performance in IC to use our ultra-pure polyethylene screw caps, see image below.



Note that these caps, are a single use cap and will not reseal after injection.

Click [HERE](#) for more information about IC compatible caps.

Click [HERE](#) for PP vials that are IC compatible.