



RSA Glass Tubes Are Not Available for Sale - Tech Information

Date: 27-APRIL-2012 Last Updated: 13-DECEMBER-2025

RSA™ glass is not available in any tube format and at this time.

RSA™ or Reduced Surface Activity glass is only available in 12x32mm autosampler vials. These vials will not release any contaminants such as chloroform or any other solvent into your sample.

For larger and smaller size vials, you may want to consider RSA-Pro X surface treatment. This is a non stick, hydrophobic surface treatment. Click [HERE](#) for information about RSA-Pro X vials and custom surface treatments.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com