

## Does MICROSOLV Sell Bared Fused Capillaries for CZE That Are Completely UV Transparent - FAQ

*Date: 11-DECEMBER-2014 Last Updated: 3-NOVEMBER-2025*

No, all of our CZE capillaries are coated with polyimide. This polymer imparts structural rigidity to the capillary but has the disadvantage of being non-UV-transparent.

Without the coating, the capillary would be very brittle and will break easily. However, this does not present a problem in typical usage since you only need a small portion of the capillary to be UV-transparent for the detection window. This can be accomplished using our [Window Maker™](#) product.

Click [HERE](#) for capillaries for CZE ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)