

Do you Offer CE Capillaries of Various Outer Diameters - FAQ

Date: 11-DECEMBER-2014 Last Updated: 3-NOVEMBER-2025

No, the only outer diameter (OD) bare fused silica capillary for ~~CE~~ that we sell is 365µm.

Other OD's may be available as CUSTOM ORDER. Minimum Order Quantities will be required. Please inquire with Customer Service. Click [HERE](#) for capillaries for ~~CE~~ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com