

What is a Zero Flow CE capillary mostly used for - FAQ

Date: 21-APRIL-2012 Last Updated: 3-NOVEMBER-2025

There are many uses for Zero Flow™ capillary in CZE. One use is to allow analytes to separate based only on electrophoretic mobility. The Zero Flow™ capillary is an excellent choice for this technique because, regardless of the pH you use to ionize your analytes, the flow will still be zero.

Another excellent use is cIEF.

Click [HERE](#) for Zero Flow and Controlled Flow capillaries ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com