

Is the Conditioning of the Controlled low and Zero Flow CE Capillaries the Same as with a Bare Fused Silica Capillary - FAQ

Date: 21-APRIL-2012 Last Updated: 3-NOVEMBER-2025

No. Inside the front cover of the capillary case is a special set of conditioning steps that must be taken.

It is important that these conditions be followed. You can also see these conditions on our website by clicking on this link: [Controlled flow capillary conditioning steps](#).

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com