



RoHS Statement for MICROSOLV CE Window Makers - Regulatory Information

Date: 12-JUNE-2015 Last Updated: 8-NOVEMBER-2025

The MICROSOLV Window Makers™ and all components used in their construction are RoHS compliant.

*Note: RoHS stands for **Restriction of Hazardous Substances Directive**. The directive mandates that RoHS compliant goods restrict the amounts of hazardous materials used in the components, namely lead, mercury, cadmium, hexavalent chromium, poly brominated biphenyls (PBB), and poly brominated diphenyl ether (PBDE). Each hazardous substance may be present in no more than 0.1% by weight for a given component (0.01% for cadmium). These limits represent the maximum thresholds for the individual components, not the final product. The directive helps ensure consumer safety by restricting the levels of these substances in electronic consumer goods.*

Click [HERE](#) for MICROSOLV Window Maker™ ordering information and pictures

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com