



Is the Power Cable Supplied with the MICROSOLV Window Maker or is it a Separate Accessory - FAQ

Date: 5-AUGUST-2015 Last Updated: 8-NOVEMBER-2025

The MICROSOLV Window Maker™ comes with a power cord that is directly built into the unit itself and is permanently hard wired. You do not need to purchase any additional power cable accessory to use this product.

Click [HERE](#) for Window Maker™ ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com