



## Does the Replaceable Heating Module Come with the MICROSOVLV Window Maker When you Order or Does it Have to be Ordered Separately - FAQ

*Date: 22-APRIL-2012 Last Updated: 1-NOVEMBER-2025*

When you order the MICROSOVLV Window Maker™, it is shipped with a heating module that creates a detection window that is approximately 2mm. The replacement module is available in the event you need to replace it or other size window modules are available in 1mm, 4mm or 7mm.

Click [HERE](#) for Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)