



Can the Window Maker Heating Modules be Used with Either US or European Voltage Units - FAQ

Date: 30-OCTOBER-2013 Last Updated: 8-NOVEMBER-2025

The Window Maker units are currently produced with 110/120V voltage for use in US & Canada (cat# [07200-S](#)) and the replacement modules (cat# [07301-S](#), [07201-S](#), [07204-S](#), [07207-S](#)) differ only in size of the burn window and can be used with new and older Window Makers even if you have a discontinued unit with different voltages.

Click [HERE](#) for MICROSOVLV Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com