



Are There Other Parts or Pieces that I Need to Use the MICROSOVL Window Maker - FAQ

Date: 22-APRIL-2012 Last Updated: 8-NOVEMBER-2025

No, the MICROSOVL Window Maker™ comes with a heating module that makes a 2mm window. No other tools are needed. Other heating modules of different window sizes can be purchased separately.

Click [HERE](#) for MICROSOVL Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com