



## What is the Voltage Requirement for the Window Maker - FAQ

*Date: 1-APRIL-2016 Last Updated: 8-NOVEMBER-2025*

**The voltage requirement for the Window Maker™ is 110/115volt at 750ma.**

Click [HERE](#) for MICROSOVLV Window Maker ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)