

## At Which Temperature Should the Zero Flow Capillaries be Stored - FAQ

*Date: 15-NOVEMBER-2013 Last Updated: 3-NOVEMBER-2025*

The capillaries should be stored at room temperature (approximately 22°C) in a dry place. The temperature these capillaries are stored at does not depend on the ID of the capillary. Many customers actually tape them to the wall near their GC instrument.



Click [HERE](#) for Zero Flow capillaries ordering information and pictures.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)