

## 20mL Vial Storage - Recommended Racks - FAQ

*Date: 4-JUNE-2013 Last Updated: 15-NOVEMBER-2025*

### **Choosing the Right Rack for 20mL Vials**

If you're storing or organizing 20mL vials, MICROSOLV offers two high-quality rack options designed for different vial types:

#### **For Standard 20mL EPA/VOA Vials**

Use catalog number [95040-SB1](#)

- Description: Storage Box/Rack with cover
- Material: Polypropylene (PP)
- Color: Violet
- Dimensions: 130mm x 130mm x 102mm
- Capacity: 16 cavities
- Fits: 20mL EPA/VOA vials

#### **For Headspace Autosampler Vials**

Use catalog number [95025-SB](#)

- Description: Storage Box/Rack with cover
- Material: Polypropylene (PP)
- Color: Blue
- Dimensions: 130mm x 130mm x 102mm
- Capacity: 25 cavities
- Fits: Headspace autosampler vials

 Click [HERE](#) for Vial Accessory ordering information and product images

Need help choosing the right rack for your application? Feel free to contact our customer support team — we're happy to assist!