

Racks for 20 ml Headspace Vial - Tech Information

Date: 4-JUNE-2013 Last Updated: 1-FEBRUARY-2026

Choosing the Right Rack for 20 mL Vials

If you're storing or organizing 20 mL vials, MICROSOLV offers two high-quality rack options designed for different vial types:

For Standard 20mL EPA/VOA Vials

Use catalog number [95040-SB1](#)

- Description: Storage Box/Rack with cover
- Material: Polypropylene (PP)
- Color: Violet
- Dimensions: 130mm x 130 mm x 102 mm
- Capacity: 16 cavities
- Fits: 20 mL EPA/VOA vials

For Headspace Autosampler Vials

Use catalog number [95025-SB](#)

- Description: Storage Box/Rack with cover
- Material: Polypropylene (PP)
- Color: Blue
- Dimensions: 130 mm x 130 mm x 102 mm
- Capacity: 25 cavities
- Fits: Headspace autosampler vials

 Click [HERE](#) for Vial Accessory ordering information and product images

Need help choosing the right rack for your application? Feel free to contact our customer support team — we're happy to assist!

AUTOSAMPLER
VIALS AND CAPS

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com