



## What Point Diameters are Available for Ultra-Fine Point Laboratory Markers - FAQ

*Date: 24-JUNE-2016 Last Updated: 8-NOVEMBER-2025*

**Three different sizes are available: 0.25mm, 0.35mm, and 0.45mm point diameter.**

The smaller diameter point will allow for finer, more detailed markings. These can be useful for situations where space to write is limited such as on a glass insert.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)