


What are the Size and Materials of the Mobile Phase Delivery Bottle Caps and O-rings - FAQ

Date: 04-APRIL-2013 Last Updated: 10-OCTOBER-2025

The mobile phase delivery bottle caps are made of polypropylene with a PTFE insert. The O-rings are made of EDPM (Ethylene Propylene Diene Monomer) and coated with FEP. The caps, inserts, and O-rings are made for use with GL 45 screw thread bottles. The ring has an outer diameter of approximately 36.4 mm and an inner diameter of approximately 28.6 mm.



Click [HERE](#) for mobile phase delivery bottle caps ordering information and pictures.