



What is the Pressure Rating of the Knitted Reactor Coils - INTERNAL ONLY

Date: 22-APRIL-2012 8-NOVEMBER-2025

The knitted reactor coils are rated to hold 300 psig. **CAUTION** : After time and extended use, this rating may be significantly reduced.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com