



What Materials are the Knitted Reactor Coils Manufactured - INTERNAL ONLY

Date: 22-APRIL-2012 Last Updated: 8-NOVEMBER-2025

The knitted reactor coils KRC are made of virgin, medical grade PTFE (a very pure form of Teflon™).

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com