

## What is the Best Way to Attach the Inlet of the Knitted Reactor Coil Tubing to the Outlet of a Tee With a PEEK Fitting that will Hold High Back Pressure - INTERNAL ONLY

*Date: 27-OCTOBER-2016 Last Updated: 8-NOVEMBER-2025*

Confirm that there is not an impediment to reasonable flow through the tubing which could increase the back pressure. We would then recommend using our items [47320-30](#) and [47320-34](#) which are high pressure PEEK HPLC tubing connectors or male nuts and ferrules.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)