

What is the Inside Diameter ID of the Opening of the Luer Lock Syringes - FAQ

Date: 21-MAY-2015 Last Updated: 3-NOVEMBER-2025

The inner diameter (ID) of the opening of the Chrom Syringes has been measured to be approximately 2.1–2.2mm.



[Syringes for Chromatography Ordering Information](#)

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com