



What is the Thickness of the SS Filter Element in the Cogent Column Filters - FAQ

Date: 22-APRIL-2012 Last Updated: 3-NOVEMBER-2025

The stainless steel filter element or "screen" is approximately 0.039" thick in the pre-column, Cogent Column Filters.

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com