



What materials are the Cogent Column Filters manufactured with - FAQ

The device (body) is made from high purity, inert PEEK with an internal stainless steel filter element and these are the only materials in the pre-column, Cogent Column Filters.

[Click here to view more information about pre-column filters for HPLC](#)



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com