



Can the Cogent Column Filters Replace my HPLC Guard Columns - FAQ

Date: 22-APRIL-2012 Last Updated: 3-NOVEMBER-2025

No, a guard column is both a particulate filter (non dissolved chemicals) and a chemical filter (removes dissolved chemicals). The [Cogent Column Filter](#) is a particulate filter and adds extra protection for your expensive analytical and guard columns.

Click [HERE](#) for column protection ordering information.



Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com