



Are the AQ Brand Syringe Filters Bi-directional in Their Flow - FAQ

Date: 22-APRIL-2012 Last Updated: 3-NOVEMBER-2025

Yes, all the AQ™ Syringe Filters have bi-directional flow capabilities except the PTFE filter as well as other syringe filters that have glass fiber pre-filters.

Although in any syringe filter format, the flow in two directions may not be equivalent because of the filter device flow design. Click [HERE](#) for AQ brand syringe filter ordering information and pictures.

Attachment:

MICROSOLV AQ filter equivalency study [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com