



Are AQ Brand Nylon Syringe Filters Compatible with DMSO - FAQ

Date: 27-NOVEMBER-2012 Last Updated: 3-NOVEMBER-2025

Yes, MICROSOVLV's AQ™ brand of syringe filters for dissolution testing and HPLC are made with a polypropylene housing which is DMSO compatible.

Our nylon membranes, both 0.45um and 0.22um are fully compatible with solutions containing DMSO. Nylon is also compatible with esters, alcohols, weak acids and bases.

Click [HERE](#) for ordering information and pictures of AQ™ brand of syringe filters.

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com