

Sample Carry Over Using One Syringe Filter for Many Samples - FAQ

Date: 22-APRIL-2012 Last Updated: 1-DECEMBER-2025

You will most likely get some kind of carry over. It is not recommended to use one filter for more than one sample.

- [Click HERE for AQ syringe filter ordering information and pictures](#)

Attachment

MICROSOLV filters equivalency study pdf [Download File](#)

Printed from the Chrom Resource Center
Copyright 2025, All Rights Apply
MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900
Fax: (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com