



For Normal Phase HPLC which syringe filter do you recommend
- FAQ

PTFE is the best membrane type for use with non polar, organic solvents required for Normal Phase Chromatography NP. Click [HERE](#) for AQ™ brand PTFE syringe filter ordering information and pictures.

Attachment

MICROSOLV filters equivalency study pdf [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com