


Is the AQ Brand Nylon Syringe Filter Compatible with Formamide - FAQ

Date: 18-February-2020 Last Updated: 3-NOVEMBER-2025


Yes, our AQ™ Brand Nylon Syringe Filters are compatible with formamide.

Formamide, dimethylformamide, and acetamide, are all compatible with our AQ™ brand Nylon Syringe Filters.


Acetamide :



Dimethylformamide:



Formamide:



Click [HERE](#) for ordering information and pictures of AQ brand of syringe filters.

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com