



## Which Solvents Should I Avoid with Nylon AQ Brand Syringe Filters - FAQ

*Date: 22-APRIL-2012 Last Updated: 3-NOVEMBER-2025*

Aldehydes and strong acids such as hydrochloric and nitric should be avoided with AQ™ Brand Nylon Syringe Filters or membrane filters.

Click [HERE](#) for syringe filters ordering information and pictures.

### Attachment

**MICROSOLV filters equivalency study pdf** [Download File](#)

Printed from the Chrom Resource Center

Copyright 2025, All Rights Apply

**MicroSolv Technology Corporation**

9158 Industrial Blvd. NE, Leland, NC 28451

Tel: (732) 380-8900

Fax: (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)